


ART FILLER UNIVERSAL LONG-LASTING PERFORMANCE PROVEN OVER 18 MONTHS

MAINTAINED CORRECTION OVER 18 MONTHS IN THE NLF

Scoring combined with profilometry



D0 before injection




D0 before injection

100% OF SUBJECTS⁽¹⁾ STILL IMPROVED



At 12 months without re-injection




At 12 months without re-injection

98,4% OF SUBJECTS⁽¹⁾ STILL IMPROVED



At 18 months




At 18 months

ART FILLER UNIVERSAL LONG-LASTING PERFORMANCE PROVEN OVER 18 MONTHS

NO HYPERECHOGENIC RIM

Evaluation by Dermal High-Frequency Ultrasound



D0 before injection

9 months

12 months with re-injection at 9 months

18 months

ART FILLER EXCLUSIVE SYRINGE⁽³⁾ FOR PRECISE AND COMFORTABLE INJECTION

OPTIMAL INJECTION EXPERIENCE



**INTEGRATED
LUER LOCK**

**HIGH VOLUME
1.2ML⁽⁴⁾**

**ERGONOMIC
FINGER GRIP
/ 360° MOBILE**

**RIGID
PLUNGER**

**WIDE
THUMB REST**

Berjaya Estetika
INDONESIA



FILLMED
LABORATOIRES

PT Berjaya Estetika Indonesia
0817-7074-2306
@fillmed_ind
@berjayaestetikaIndonesia

www.berjayaestetika.id

ART FILLER HAUTE-SCULPTURE CORRECTION

TRIAL-HYAL INNOVATION

3 HA - NATURAL FLEXIBLE ENTANGLING
3 hyaluronic acids, 3 correction potentials
(sculpting, volumising, smoothing), optimal mouldability

EASY-TO-SHAPE TEXTURES

MATCHING EACH INJECTION SITE
Optimised texture tailored to match each injection site:
superficial lines, lips, medium to deep wrinkles, volume restoration

HIGH PERFORMANCE & TOLERANCE

PROVEN OVER 18 MONTHS⁽¹⁾
Correction maintained and high tolerance was proven during
the clinical evaluation.
Bling evaluation and instrument - based measurements⁽²⁾
(profilometry and dermal high-frequency ultrasound imaging)

OPTIMAL COMFORT

FOR OPTIMAL INJECTION EXPERIENCE
Products with lidocaine 0.3% high injectability.
Ergonomic syringe

MADE IN EUROPE: FILLMED MANUFACTURING MAXIMAL QUALITY

International standards ISO 14 644, ISO 13 485

TRY-HYAL

	WRINKLES		VOLUME	LIPS	
	FINE LINES	UNIVERSAL	VOLUME	LIPS SOFT	LIPS
INDICATION	SOFT & SMOTHING Fine and superficial wrinkles	FILLING & SCULPTING Medium to deep wrinkles	VOLUME & SHAPING Volume creation & restoration	SCULPTING & SMOOTHING Perioral wrinkles and soft volume	PLUMPING & DEFINING Lip volume
HA CONCENTRATION	20 mg/mL	25 mg/mL	25 mg/mL	25 mg/mL	25 mg/mL
RATIO: <small>VERY LONG CHAINS</small> LONG CHAINS	●○○○	●●○○	●●●●	●●○○	●●●○
CROSS-LINKING DEGREE	●○○○	●●●○	●●●●	●●●○	●●●○
FREE HA	●●●●	●●●○	●○○○	●●●○	●●○○
LIDOCAINE	0.3%	0.3%	0.3%	0.3%	0.3%
PACKAGING	2 x 1 mL	2 x 1.2 mL	2 x 1.2 mL	2 x 1.2 mL	2 x 1.2 mL
NEEDLES	4 x 30G½"	4 x 27G½"	4 x 27G½"	4 x 30G½"	4 x 27G½"




ART FILLER TRI-HAL INNOVATION NATURAL FLEXIBLE ENTANGLING

3 HYALURONIC ACIDS

- 1-VERY LONG CHAINS
- +
- 2-LONG CHAINS



Network of
natural flexible entangling
> sculpting & volumising
potentials




3-FREE HYALURONIC ACID
> smoothing potential

3 CORRECTION POTENTIALS



1-SCULPTING POTENTIAL
EASY-TO-SHAPE: capacity
of the gel to be sculpted



2-VOLUMISING POTENTIAL
ELASTICITY: capacity of
the gel to resist to stress



3-SMOOTHING POTENTIAL
VISCOSITY: capacity of the
gel to integrate smoothly
into the tissues